





ALL SHOPS

304 STAINLESS STEEL




STATIC DISSIPATIVE PVC WINDOW

WALL CUT OUT 12.500W X 12.500H



ELECTROPOLISHED 304 STAINLESS STEEL CHAMBER

REVISIONS				MATERIAL SPEC.	FINISH SPEC.	UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES	TITLES
RV#	DESCRIPTION	DATE	BY				
NC	INITIAL RELEASE	4-26-18	JR	STAINLESS STEEL		TOLERANCES ON: .XX DECIMALS ± .03 X" ± 1" .XXX DECIMALS ± .010 X/X ± 1/64	 Terra Universal PASS THROUGH CHAMBER 304 SS WALL MOUNT 12Wx12Dx12H
○				SD PVC			

Proprietary & Confidential
All rights reserved by
TERRA UNIVERSAL
Disclosure or Copying Prohibited