





ALL SHOPS

WALL CUT OUT
24.500W X 24.500H


SECTION U-U



mount bracket loose





SPRING LOADED
SS HINGE



UL RATED
FIRE RATED GLASS

SS PADDLE
LATCH
AS07478


FIRE RATED
DOOR



SS LIFT LATCH

STATIC PVC
WINDOW

SS HINGE (X2)

REVISIONS				MATERIAL SPEC.	FINISH SPEC.	UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES	 <p>Proprietary & Confidential All rights reserved by TERRA UNIVERSAL Disclosure or Copying Prohibited</p>	TITLE	
RV#	DESCRIPTION	DATE	BY	304SS		TOLERANCES ON:		<p>PASS-THROUGH, 304 SS FIRE RATED DOOR 24"W X 24"H X 24"D</p>	
NC	INITIAL RELEASE	02-18-14	JR			.XX DECIMALS ± .03 X" ± 1"			
①	SHEET 28 LIFT LATCH PHOTO	8-24-16	G.F.			.XXX DECIMALS ± .010 X/X ± 1/64			
	SHEET 29 INTERLOCK					SIGNATURES	DATE		
				ISSUED BY: PHILLIP PAVLOVICH	6/10/2019			A	2641-12
				REF #: 195844				SCALE	SHEET 1 OF 30