

# ADA-Compliant Protector Solvent Storage Cabinet with self-closing door(s)

<http://www.labconco.com/product/ada-compliant-protector-solvent-storage-cabinet-with-self-closing-doors/1040>

Catalog Number: 9906200



**Weight:** 120.0 lbs

**Weight metric:** 54.0 kg

**Dimensions:** 24.0" w x 22.0" d

**Dimensions metric:** 61.0 x 55.9 cm

**Cabinet Width:** 24" (61 cm)

**Door(s):** Single, Right Hinge, Automatic Close

**Filler Panel Depth:** 8.0" (20.3 cm)

**Maximum Base Height:** 32.75" (83.2 cm)

**Minimum Base Height:** 31.5" (80.0 cm)

**Storage Type:** Solvent storage

Protector Solvent Storage Cabinets support Labconco fume hoods and safely store and vent solvents and other flammable liquids. In case of a building fire, these cabinets protect the contents inside them.

## Features

- Support Labconco fume hoods including Protector Premier, XStream, ClassMate, Stainless Steel, and PVC; Protector XL Benchtop and HOPEC IV; Basic and Fiberglass 30 Hoods
- Durable dual wall epoxy-coated steel construction with 1.5" air gap
- Attractive glacier white exterior complements laboratory casework and Labconco fume hoods.
- Include epoxy-coated steel shelf
- Include two rear vent connections with flame arrestors and closure plugs
- Include four leveling feet and one 8" filler panel to increase cabinet depth from 22" to 30".
- Support loads up to 800 pounds.

## Options

- Manual-closing door(s) with 3-position locking handle
- Self-closing door(s) with 3-position locking handle. Self-closing door(s) close automatically in case of fire when temperature reaches 74° C (165° F)

## Required Accessories

- SpillStopper™ Work Surfaces
- Flat Work Surfaces

## Optional Accessories

- Shelf Kits
- Filler Panels

Continue reading online at <http://www.labconco.com/>

## For Product Assistance



**Dan Simpson**  
Sales Representative  
1-800-821-6699,  
extension 109



**Steve Marchand**  
Senior Sales  
Representative  
1-800-821-6699,  
extension 111



Labconco Corporation  
8811 Prospect Avenue  
Kansas City, MO 64132-2696  
Phone: 816-333-8811  
Toll Free: 800-821-5525

Scan for more info:

