

Information System Component Mounting Solutions Flat Panel/Keyboard Combinations

Ergo-Bracket



L-Bracket



Folding Keyboard Ergo-Bracket



Provides recommended ergonomic separation between the display and keyboard.
Provides tilt adjustment for the flat panel for optimal viewing angle.

Example Configurations

Keyboard Folding Ergo-Bracket



Ergo-Bracket



L-Bracket



Ergo Bracket for Flat Panel Keyboard Combinations FLP-0004-83

- Provides a secure, ergonomic mounting solution for both the flat panel display and keyboard
- Provides ergonomic separation of the display and keyboard with independent display tilt adjustment
- Adjustable tilt for display
- VESA® FDMI™ 75/100 mm compatible
- Accepts all GCX keyboard trays/plates
- Product weight: 7.5 lbs/3.4 kg

L-Bracket for Flat Panel/Keyboard Combinations FLP-0004-84

- Provides a compact mounting solution for the flat panel display and keyboard
- Flat panel display and keyboard mounts
- Display height is adjustable to accommodate larger displays and user preferences
- Tilt adjustable (display and keyboard in unison)
- VESA® FDMI™ 75/100 mm compatible

Folding Ergo-Bracket for Flat Panel/Keyboard Combinations FLP-0004-89

- Folding ergo-bracket allows the keyboard tray/plate to be folded upward for security or to save space
- Provides ergonomic separation of display and keyboard with independent display tilt adjustment
- Accepts all GCX keyboard trays/plates
- VESA® FDMI™ 75/100 mm compatible

Keyboard Mounts

Keyboard Tray – 20.8"/52.8 cm with Left/Right Mouse Trays WM-0023-47

- 20.8"/52.8 cm wide tray for standard keyboard with left/right 8"/20.3 cm x 6.5"/16.5 cm slide-out mouse trays
- Padded wrist rest and side bumpers
- Includes mouse house, mouse pads

Keyboard Tray – 25.8"/65.5 cm WM-0023-32

- 25.8"/65.5 cm wide tray
- Padded wrist rest and side bumpers
- 8"/20.3 cm x 8"/20.3 cm mouse surfaces (both sides)
- Includes mouse house and mouse pads
- Product weight: 5.5 lbs/2.6 kg



GCX is an ISO 9001 registered company (TUV-Essen)
GCX, VHRC, VHRS, VHM, PRC, M-Series, PolyMount and PolyQuip are trademarks of GCX Corporation.

© Copyright 2010

MKT-0001-08 2/2010