



## PROHT

- All stainless steel construction
- 14 gauge 316L stainless steel tank, 304 stainless steel skirt
- Powered by Neptune® C 40kHz generator with sweep frequency
- Digital controls for the heat and ultrasonics
- Ultrasonic frequency sweep to enhance cleaning
- Uniform cavitation with rapid degassing
- Basket
- Lid
- Voltage: 120V or 240v, single phase



## PROHT - Series

1206 / 2406 / 3606 / 6006 / 1212 / 2412 / 3612  
1812 / 2016 / 2020 / 2424 / 3624



CLEANING TECHNOLOGIES GROUP

RANSOHOFF • BLACKSTONE-NEY ULTRASONICS • CTG ASIA

**Blackstone-NEY Ultrasonics**  
(800) 766-6606  
P.O. Box 220  
9 N. Main St.  
Jamestown, NY 14702-0220  
[www.ctgclean.com](http://www.ctgclean.com)

**Ransohoff**  
(800) 248-9274  
4933 Provident Drive  
Cincinnati, OH 45246

**CTG Asia**  
(86)0512-66161698  
56 Songshan Road  
Suzhou New District  
Jiangsu Province, China PRC: 215151



# CLEANING TECHNOLOGIES GROUP

RANSOHOFF • BLACKSTONE-NEY ULTRASONICS CTG ASIA

MODEL NUMBER	PROHT-1206	PROHT-2406	PROHT-3606	PROHT-6006	PROHT-1212	PROHT-2412
VOLTAGE	120V	120V	120V	120V	240V	240V
CAPACITY	1.88 Gal (7.1 Ltrs)	3.7 Gal (14 Ltrs)	5.6 Gal (21.2 Ltrs)	9.4 Gal (35.6 Ltrs)	5.6 Gal (21.2 Ltrs)	11.2 Gal (42.4 Ltrs)
TANK DIMENSIONS (LxWxH)	305mm x 152mm x 203mm	610mm x 152mm x 203mm	914mm x 152mm x 203mm	1524mm x 152mm x 203mm	305mm x 305mm x 279mm	610mm x 305mm x 279mm
WORKING ENVELOPE (LxWxH)	12"x6"x8"	24"x6"x8"	36"x6"x8"	60"x6"x8"	12"x12"x11"	24"x12"x11"
UNIT DIMENSIONS (LxWxH)	457mm x 305mm x 356mm	762mm x 305mm x 356mm	1067mm x 305mm x 356mm	1676mm x 305mm x 356mm	457mm x 457mm x 432mm	762mm x 457mm x 432mm
UNIT DIMENSIONS (LxWxH)	18"x12"x14"	30"x12"x14"	42"x12"x14"	66"x12"x14"	18"x18"x17"	30"x18"x17"
ULTRASONIC POWER	240 Watts	600 Watts	960 Watts	1,440 Watts	480 Watts	960 Watts
HEAT POWER	600 Watts	1,200 Watts	2,000 Watts	3,000 Watts	1,200 Watts	2,000 Watts

MODEL NUMBER	PROHT-3612	PROHT-1812	PROHT-2016	PROHT-2020	PROHT-2424	PROHT-3624
VOLTAGE	240V	240V	240V	240V	240V	240V
CAPACITY	16.8 Gal (63.6 Ltrs)	11.2 Gal (42.4 Ltrs)	16.6 Gal (62.9 Ltrs)	24.2 Gal (91.6 Ltrs)	34.9 Gal (132.1 Ltrs)	52.4 Gal (198.4 Ltrs)
TANK DIMENSIONS (LxWxH)	914mm x 305mm x 279mm	457mm x 305mm x 356mm	508mm x 406mm x 356mm	508mm x 508mm x 406mm	610mm x 610mm x 406mm	914mm x 610mm x 406mm
WORKING ENVELOPE (LxWxH)	36"x12"x11"	18"x12"x14"	20"x16"x14"	20"x20"x16"	24"x24"x16"	36"x24"x16"
UNIT DIMENSIONS (LxWxH)	1067mm x 457mm x 432mm	610mm x 457mm x 508mm	660mm x 559mm x 508mm	660mm x 660mm x 559mm	762mm x 762mm x 559mm	1067mm x 762mm x 559mm
UNIT DIMENSIONS (LxWxH)	42"x18"x17"	24"x18"x20"	26"x22"x20"	26"x26"x22"	30"x30"x22"	42"x30"x22"
ULTRASONIC POWER	1,680 Watts	960 Watts	960 Watts	1,920 Watts	2,880 Watts	3,840 Watts
HEAT POWER	3,000 Watts	2,000 Watts	3,000 Watts	4,000 Watts	5,000 Watts	7,000 Watts